

Sertraline

Description	Zoloft®, is used for the symptomatic relief of depressive illness. However, the antidepressant action of sertraline in hospitalized depressed patients has not been adequately studied.
Therapeutic Low	0.0100
Therapeutic High	0.2000
Units	ug/mL
Specimen	P
Half Life Low	22.00
Half Life High	36.00
VOD Low	76.00
VOD High	0.00
Toxic Level	1.200
Lethal Level	1.500
Plasma to Whole Blood ratio	0.00
Warnings	Warnings - may impair mental and/or physical ability required for the performance of potentially hazardous tasks (e.g., driving, operating heavy machinery).

Desmethylsertraline

Description	Desmethylsertraline is the predominant metabolite of the antidepressant sertraline, Zoloft®. While it is an active metabolite, it is substantially less active than sertraline.
Therapeutic Low	0.0000
Therapeutic High	0.0000
Units	ug/mL
Specimen	
Half Life Low	62.00
Half Life High	104.00
VOD Low	0.00
VOD High	0.00
Toxic Level	0.000
Lethal Level	0.000
Plasma to Whole Blood ratio	0.00
Warnings	Warnings - may impair mental and/or physical ability required for the performance of potentially hazardous tasks (e.g., driving, operating heavy machinery).